

Processing techniques

Depending on the processing technique and indication, the following recommendation can be given with regard to the selection of materials.

	FULL CONTOUR					CUT-BACK				LAYERING		
	Inlay / Onlay	Veneer	Partial Crown	Crown	Bridge	Veneer	Partial Crown	Crown	Bridge	Veneerable inlay	Coping	Bridge
RESIN NANO CERAMIC												
3M™ ESPE™ Lava™ Ultimate Restorative – LT & HT	•	•	•	•								
CERAMIC												
3M™ ESPE™ Lava™ Plus Zirconia ¹				•	•			•	•	•	•	•
zerion®										•	•	•
IPS e.max® CAD – HT	•	•	•	•		•	•	•				
IPS e.max® CAD – LT			•	•		•	•	•				
IPS e.max® CAD – MO											• ²	
IPS Empress® CAD – HT	•	•				•						
IPS Empress® CAD – LT			•	•	•	•		• ¹				
IPS Empress® CAD – Multi			•	•	•	•						
VITA Mark II	•	•	•	•		•		• ¹				
VITA TriLuxe			•	•	•	•		• ¹				
METAL												
ticon®				•	•			•	•	•	•	•
coron®				•	•			•	•	•	•	•
POLYMERS												
polycon® ae				•	•			•	•	•	•	•
polycon® cast				•	•			•	•	•	•	•

• Recommended

• Possible

¹ Only for anterior crowns

² Only for anterior crowns and premolars