

Cementing

Depending on the indication and selected material, the CARES® tooth restorations can be fitted using adhesive, self-adhesive or conventional cementing (see table below).

	Adhesive	Self-adhesive	Conventional	Temporary
RESIN NANO CERAMIC				
3M™ ESPE™ Lava™ Ultimate Restorative – LT & HT	●	●		
CERAMIC				
3M™ ESPE™ Lava™ Plus Zirconia ¹	●	●	●	
zerion®	●	●	●	
IPS e.max® CAD	●	● ¹	● ¹	
IPS Empress® CAD	●			
VITA Mark II & Triluxe ¹	●	● ¹		
METAL				
ticon®			●	
coron®			●	
POLYMERS				
polycon® ae				●
polycon® cast	n. a.			

¹ For crowns only