

# **SINGLE-PATIENT DRILLS**



COMMITTED TO  
**SIMPLY DOING MORE**  
FOR DENTAL PROFESSIONALS

## FEATURES AND BENEFITS

### Reliable

- Drill sharpness and optimal cutting properties are ensured with single-patient drills.
- The drills are delivered in a sterile package and offer improved protection of patients and practice personnel against the risk of infection.

### Simple

- Time-saving as single-patient drills do not require cleaning or sterilization prior to use.
- Clear diameter color-coding on the drills facilitates the surgical procedure.
- Single-patient drills have the same dimensions, cutting properties and depth marks as Straumann multi-use drills.
- Provides for easier cost management per case.
- Easy handling due to integrated sterile carrier unit that allows aseptic delivery of the drills onto sterile field by a non-sterile person.

### Versatile











- For the different implant diameters, three different drill sets are available for order. Drills with a diameter of 3,5 mm or 4,2 mm can be ordered individually.
- All drills are available in a short and a long version.



Single-patient drills represent a safe and simplified treatment solution, especially with regard to patient and case management. They facilitate cost assessment since each drill is only used once per patient and surgery.

# OPTIMAL CUTTING PROPERTIES AND SIMPLIFIED TREATMENT COST MANAGEMENT

## Overview of the different drill sets

Art. No.		Article	Dimensions	Material
040.400S		Drill set, short, for Ø 3,3 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, length 33 mm	stainless steel
040.401S		Drill set, long, for Ø 3,3 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, length 41 mm	stainless steel
040.403S		Drill set, short, for Ø 4,1 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, Ø 3,5 mm, length 33 mm	stainless steel
040.404S		Drill set, long, for Ø 4,1 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, Ø 3,5 mm, length 41 mm	stainless steel
040.406S		Drill set, short, for Ø 4,8 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, Ø 3,5 mm, Ø 4,2 mm, length 33 mm	stainless steel
040.407S		Drill set, long, for Ø 4,8 mm implants, single-patient use, sterile	Ø 2,2 mm, Ø 2,8 mm, Ø 3,5 mm, Ø 4,2 mm, length 41 mm	stainless steel
044.234S		Twist drill, short, for Ø 4,1 mm implants, single-patient use, sterile, red	Ø 3,5 mm, length 33 mm	stainless steel
044.235S		Twist drill, long, for Ø 4,1 mm implants, single-patient use, sterile, red	Ø 3,5 mm, length 41 mm	stainless steel
044.238S		Twist drill, short, for Ø 4,8 mm implants, single-patient use, sterile, green	Ø 4,2 mm, length 33 mm	stainless steel
044.239S		Twist drill, long, for Ø 4,8 mm implants, single-patient use, sterile, green	Ø 4,2 mm, length 41 mm	stainless steel

### Indication

Single-patient drills are indicated for all implant bed preparation procedures carried out with the Straumann® Dental Implant System. They are supplied in a sterile package and must not be used on more than one patient and in more than one surgery. Single-patient drills can help minimize the risk of patient infection and are color-coded for simplified identification.

[www.straumann.com](http://www.straumann.com)



STRAUMANN GUARANTEE

**International Headquarters**

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01